

SAFETY BOX WITH FOAM INSERT FOR KLESSIDRA (21 places)



CODE	DESCRIPTION
05-900700	Safety box complete with foam with 21 places – 1 pcs
05-900800	Accessory: foam with 21 places – 1 pcs

Distributor: Bio-Optica Milano S.p.A.

Date of emission: 13/02/2020

Rev. 00

DESCRIPTION

Safety box for Klessidra 1.0 and 2.0 complete with foam support with 21 places.

Made of polycarbonate. The cover has a hermetic silicon seal, to ensure a high degree of safety, even in case of accidental leak of liquids contained within, and it is provided with 4 safety locks preventing opening if the box drops down.

Complete with stainless steel handle. Biohazard logo printed on the lid.

TECHNICAL FEATURES

Safety box	Dimensions	Height	180 mm
		Length	330 mm
		Depth	175 mm
	Material	Polycarbonate, colourless	
Foam support	Dimensions	Height	55 mm
		Length	330 mm
		Depth	140 mm
		Hole diameter	35 mm
		Place number	21
	Material	Moltoprene, dark grey	

Date of emission: 13/02/2020

Rev. 00