

Accessories for archives



CODE	DESCRIPTION	PACKAGING
03-5000-MDL	Separating spring	4 pcs
03-C28N	Basin for slides and blocks	1 pcs
03-C28S	Basin for mega slides and supermegacassettes	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of emission: 17/03/2020

Rev. 02


Accessories for big archiving systems:

- Separating spring for slides;
- Basin for slides and blocks.

SLIDES SPRING FEATURES



The metallic separating spring is intentionally realized to separate slides still wet (just mounted), avoiding that they remain attached one to the others. Each spring is long 15 cm, with a diameter of 2,5 cm and it is perfectly contained in a compartment of C28 basin. Each spring may contain until a maximum of 48 slides.



BASIN FEATURES



The basins are realized in plastic, grey coloured, to contain slides and blocks in the big archives. They are available in three different models:

- Basin C28N: standard basin. External dimensions 30 x 6,35 x 4 (h) cm.

Basins model C28N are divided in four compartments, each of these may contain until 20 cassettes (total 80 cassettes for each basin) or 120 slides (total 480 slides for each basin). In each compartment, a spring may be inserted: in this case, the capacity is reduced at 48 slides (total 192 slides for each basin equipped with 4 springs).

- Basin C28S: specifically studied for archiving supermegacassettes and mega-slides, without internal compartments. The external dimensions are the same of basin C28N (30 x 6,2 x 4 (h) cm). Each basin C28S may contain 9/10 supermegacassettes or 240 mega-slides.

Date of emission: 17/03/2020

Rev. 02

Technical specifications

Spring	Material	Stainless steel
	Dimensions (l x Ø)	(150 x 25) mm
Basin	Material	Polystyrene
	Dimensions (l x d x h)	(300 x 63,5 x 40) mm
Packaging	Primary packaging	Cellophane bag
	Label	Adhesive label on external packaging
	Secondary packaging	N.a.
Storage	Storing conditions	Not expected
	Storage temperature	Not expected
	Stability	The product is stable over the time
	Validity	Not applicable
Warning and precautions	Mode of use	Precautions not expected
	Classification	For specialized technical operators, not classified as dangerous product
	Disposal	Dispose in conformity with local laws

Date of emission: 17/03/2020

Rev. 02