




**MICROSCOPE SLIDES**  
*for big sections*



CODE	BAND	EDGE	CORNER	PACKAGING
09-7000	Double frosted	Ground edge 90°	Clipped 45°	1.250 pcs (25 boxes x 50 slides)



In vitro diagnostic – medical device



Manufacturer: Bio-Optica Milano S.p.A.



Disposable product

Date of issue: 19/07/2021  
Rev. 01

### TECHNICAL FEATURES

Microscope slides for big sections.  
 Dimensions (52) mm x (76) mm  
 Thickness: 1-1.1 mm  
 Height band: 19.5(±0.5) mm

Composition:	Extra white glass <b>in accordance with ISO 8037 standard</b> <ul style="list-style-type: none"> <li>• SiO<sub>2</sub>: 72.30%</li> <li>• Na<sub>2</sub>O: 14.20%</li> <li>• K<sub>2</sub>O: 1.20%</li> <li>• CaO: 6.40%</li> <li>• MgO: 4.30%</li> <li>• Al<sub>2</sub>O<sub>3</sub>: 1.20%</li> <li>• Fe<sub>2</sub>O<sub>3</sub>: 0.03%</li> <li>• SO<sub>3</sub>: 0.30%</li> </ul>
Refraction index:	1.513 – 1.523 (measured between $\lambda= 546.07$ nm and $\lambda= 643.85$ nm)
Density:	(2.47 ± 0.01) kg/dm <sup>3</sup>
Primary packaging:	50 carton boxes with 25 slides each
Secondary packaging:	Carton box
Validity:	3 years
Storage:	Keep the slides in a fresh and dry environment. Avoid large variations in temperature during both storage and usage. Cooling of the product can cause condensation and lead to condensed water forming between the glasses
Disposal:	Dispose in conformity with the local law. If not used, dispose like common glass.

Date of issue: 19/07/2021  
 Rev. 01