

## **STAINLESS STEEL EMBEDDING MOLDS**



CODE	DIMENSIONS	PACKAGING
07-BM15155	15 x 15 x 5 mm	12 pcs
07-BM24245	24 x 24 x 5 mm	12 pcs
07-BM30245	30 x 24 x 5 mm	12 pcs
07-BM37245	37 x 24 x 5 mm	12 pcs
07-BM775	7 x 7 x 5 mm	12 pcs
07-MBM6	33 x 24 x 12 mm	6 pcs
07-7010	65 x 45 x 15 mm	10 pcs



In Vitro Diagnostic – Medical Device



Manufacturer: Bio-Optica Milano S.p.A.

Date of emission: 20/03/2020

Rev. 03

Varied measures steel tiles for inclusion of histological paraffin samples.

### TECHNICAL SPECIFICATIONS

Bio Molds are reusable products. To guarantee the cleaning, we advise to clean them with xylene or d-limonene, after use, to eliminate all paraffin remnants.

Technical data	Composition	Stainless Steel
	Dimensions	See table at page 1
Packaging	Primary container	Cellophane bag
	Secondary container	None
Storage	Storing conditions	None
	Validity	Not applicable to this kind of product.
Warning and precautions	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of emission: 20/03/2020

Rev. 03