

## ***Cartoglass and Cartobloc*** Cardboard archives



CODE	DESCRIPTION	PACKAGING
03-4010-BC	Cartobloc, cardboard archive for blocks	10 pcs
03-4015-BA	Cartoglass, cardboard archive for glass slides	10 pcs
03-4020-CO	Cover for Cartoglass and Cartobloc	10 pcs



Manufacturer: Bio-Optica Milano S.p.A.  
Via San Faustino 58  
20134 MILANO (MI)

Date of issue: 22/10/2018  
Rev. 02

Archives for slides and histological blocks, realized in resistant, but light, cardboard. They are easy to transport, also if they are full.

The structure is realized in light and resistant cardboard, with internal separators. The slides archive Cartoglass can contain up to 3000 slides; the block archive Cartobloc can contain until 320 blocks or 220 rings.

### Technical Specifications

Technical data	Material	Cardboard
	Dimensions Cartobloc (l x p x h)	290 x 400 x 45 mm
	Dimensions Cartoglass (l x p x h)	290 x 400 x 80 mm
	Dimensions Cover (l x p x h)	290 x 400 x 15 mm
	Capacity Cartobloc (blocks)	320
	Capacity Cartoglass (slides)	3000
Packaging	Primary container	Cardboard box
	Lable	Paper adhesive lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 22/10/2018

Rev. 02