

## Slides mailer

### General features

Slides mailers.



Code	Description	Dimension	Packaging
09-000530	For 5 slides, with pressure cap	29 x 82 x 16 mm	50 pcs.
44-13061-50	For 5/10 slides, with screw cap	∅ 40 x 90 mm	50 pcs.
44-13031	For 1 slide	50 x 100 x 6 mm	500 pcs.
44-13041	For 2 slides	73 x 85 x 6 mm	500 pcs.
44-13051	For 3 slides	100 x 84 x 6 mm	100 pcs.

Distributed by Bio-Optica Milano S.p.A.

### Technical features

Raw material: Polypropylene.

Temperature range: -10°C / + 120°C.

These slides mailers are ideal for the dispatching or storage of small quantities of slides. Thanks to the thickness of the walls, slides are fully protected from accidental breakage. In the slides mailer with screw cap (code 44-13061-50), the slide protrudes about 20 mm from the superior board: so, it could be easily extracted.