

Cover for Archives



In the photo, from the top:

- Cover
- Module for macrosections
- Module Colorteca
- Module Istobloc
- Module Istoglass
- Base

| CODE | DESCRIPTION | PACKAGING |
|------------|--|-----------|
| 03-5000-CO | Cover for modular archives Istoglass, Istobloc and Colorteca | 1 pcs |



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 22/10/2018
Rev. 02

Cover for modular archives, realized in sheet metal, thickness 9/10, japanned white.
It has to be used with:

- Istoglass, code 03-5000-14;
- Istobloc, code 03-B900;
- Colorteca, series code 03-COLOR13.

Technical Specifications

| | | |
|-------------------------|-------------------------------|--|
| Technical data | Material | Metal sheet |
| | Dimensions (l x p x h) | 490 x 490 x 25 mm |
| | Electrical wires | Not necessities |
| Packaging | Primary container | Cardboard box |
| | Lable | Paper adhesive lable on the external packaging. |
| | Secondary container | Not used. |
| Storage | Storing conditions | No particular conditions needed. |
| | Temperature | No particular conditions needed. |
| | Expiry date | Not applicable for this kind of product. |
| | Validity | Not applicabile. |
| Warning and precautions | Mode of use | No precautions needed for the use of the product. |
| | Classification of the product | For specialized technical operators. Product not classified as dangerous. |
| | Disposal | Dispose in conformity with local laws. |

Date of issue: 22/10/2018

Rev. 02