

## **Macrosections Istoteca**

Archive for slides and blocks of big dimensions



*In the photo, from the top:*

- Cover
- Module for macrosections
- Module Colorteca
- Module Istobloc
- Module Istoglass
- Base

CODE	DESCRIPTION	PACKAGING
03-7000	Archive for macrosections (microscope slides and blocks)	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.  
Via San Faustino 58  
20134 MILANO (MI)

Date of issue: 22/10/2018

Rev. 02

The filing system Istoteca for macro sections is an archive realized in sheet metal, thickness 9/10, japanned white. Each Istoteca has seven drawers running on guides. Each drawer is furnished with a knob and the possibility to apply a paper label.

Each filing system may contain:

- about 220 Super mega cassettes for macro-sections (dimension about 67 x 76 x 15 mm);
- about 5000 slides for macro-sections (slide maximum dimension 67 x 76 mm);
- about 1645 film slides (dimension max. 67 x 70 x 2 mm), 235 for each drawer.

The size is the same of standard Istoglass and Istobloc. In this way it is possible to overlap them on standard filing systems. Each filing system can be equipped with a base (code 03-5000-BA) and a cover (03-5000-CO). They are made of metal sheet, thickness 9/10, japanned white.

It is possible overlap different Istoglass elements (until maximum 10 elements). In this way, it is possible decreasing the archives' size, assuring the content protection.

### Optional accessories

CODE	DESCRIPTION	PACKAGING
03-5000-BA	Base for Istoglass, white colour	1 pcs
03-5000-CO	Cover for Istoglass, white colour	1 pcs
03-5000-MD	Metal spring	1 pcs

### Technical Specifications

Technical data	Material	Metal sheet
	Dimensions (l x p x h)	490 x 490 x 140 mm
	Weight (empty)	13 kg
	Weight (full)	42 kg
	Drawers per module	7
Packaging	Primary container	Cardboard box
	Lable	Paper adhesive lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 22/10/2018

Rev. 02