

Separating spring for Istoglass



CODE	DESCRIPTION	PACKAGING
03-5000-MD	Separating spring for modular archive for microscope slides Istoglass	1 pcs



Manufacturer: Bio-Optica Milano S.p.A.
Via San Faustino 58
20134 MILANO (MI)

Date of issue: 22/10/2018
Rev. 02

Metal spring to use with modular archives for slides Istoglass (code 03-5000-14.)

It is ideal to separate slides when they are not completely dry, avoiding that they remain attached one to another.

In each drawer with one spring it is possible to archive up to 145 slides.



Technical Specification

Technical data	Material	Stainless
	Length	430 mm
	Diametre	25 mm
	Capacity	145 slides
Packaging	Primary container	Polyethylene bag
	Lable	Paper adhesive lable on the external packaging.
	Secondary container	Not used.
Storage	Storing conditions	No particular conditions needed.
	Temperature	No particular conditions needed.
	Expiry date	Not applicable for this kind of product.
	Validity	Not applicabile.
Warning and precautions	Mode of use	No precautions needed for the use of the product.
	Classification of the product	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of issue: 22/10/2018

Rev. 02