

## **RACK FOR KLESSIDRA (16 places)**



CODE	PACKAGING
05-900900	Shockproof plexiglass rack for Klessidra – 2 pcs



Manufacturer: Bio-Optica Milano S.p.A.

Date of issue: 28/06/2018  
Rev. 01

### PRODUCT INFORMATION

Transparent plexiglass rack with 16 places for an easy transport of Klessidra.  
The rack is provided with two handles for the transport.

### TECHNICAL FEATURES

Dimensions	Height	65 mm
	Length	223 mm
	Depth	199 mm
	Holes diameter	35,5 mm
	Plexiglass thickness	4 mm
	Place number	16
Material	Shockproof transparent plexiglass	
Storage	Temperature	Not applicable
	Expiry date	Not applicable

Date of issue: 28/06/2018

Rev. 01