

Microscope Slides Boxes

General features

In shock-resistant polystyrene, with numbered slots for 26 x 76 mm slides. They are supplied with an index-card for the classification of the prepared slides. Lids are totally removable. The 50 and 100 place boxes have locating lugs so boxes can be stacked.

Virtually unbreakable.

Code	Number of slides	Dimension (mm)	Packaging
44-13071	25 slides	98x 83 x 38 (h)	1 pcs
44-13072	50 slides	230x 97 x 38 (h)	1 pcs
44-13073	100 slides	230x 180 x 38 (h)	1 pcs



Distributed by Bio-Optica Milano S.p.A.

Technical features

- * Microbiological status: not sterile;
- * Raw material: shock-resistant polystyrene;
- * Temperature range: min -10°C, max +70°C;
- * Colour: grey;
- * Slides dimensions (mm): 26 x 76;
- * Shelf life: 5 years

The product doesn't request CE mark. The product is only for professional use in analysis laboratories.